

Our mission is to improve the quality of life for those who are visually impaired through awareness, resources and solutions

Dr. Lynne Noon
Diplomate in Low Vision Rehabilitation

Dr. Corona Hoang

The Traveller

Optelec's Revolutionary New Video Magnifier

The **Traveller HD** allows you to comfortably read letters, newspapers and other material at the dining table, in bed or in your favorite chair. Wherever you are, at home, in the office, at school or on vacation, the **Traveller HD** is your ideal light weight portable reading companion.

Improve your reading fluency with the **Traveller HD**. Its unique *Slide & Read* mechanism provides you with a more natural reading experience, especially for those people that find it difficult to keep their hands steady. Simply slide the screen across a line of text and view the content on-screen. It has never been easier to follow a line or sentence. Sliding back and forth enables you to recognize the shapes of individual words. You can see a more complete line of text because of the large screen, so you can read even more fluently. Place the **Traveller HD** on the document you want to read. To navigate vertically, use the roller and move the unit up and down over the document. To read or to navigate horizontally, slide the screen and move it from left to right, and from right to left to move back to the beginning of the line. You can work with the Optelec **Traveller** at either ViewFinder Low Vision Resource Center or call for an in-home demonstration.



The large 13.3 inch screen makes it easier to recognize the words that are being displayed. The brightness can be customized from ultra bright to soft and sharp.



The roll enables you to effortlessly navigate from the top to the bottom of the document

The slides help you to follow a straight line.

The Zoom wheel to adjust the level of magnification up to 30X

Press the Open / Close button to unlock the Traveller HD and lock it securely into place. Carrying case included for easy travel.

New BRIGHT Head Light LED

Superior Quality, Versatile Design, Attaches Easily to Eyewear!



Eschenbach has developed a compact LED lighting system that weighs just 2 ounces and will fit onto spectacles, telescopic systems, and magnifiers. Two bright LEDs provide an outstanding hands-free illumination for reading, crafts and other near tasks.

The **Head Light LED** uses 3 AA batteries (included) and the operating time is about 8 consecutive hours in full LED bright mode. Using the Eco mode (about 75% brightness), the operating time is approximately 13 consecutive hours.

The batteries are contained in a sleek case that can be placed in a pocket or on a belt using the optional belt clip. The battery case is connected to the LEDs by a thin 3.2 foot cable that can be attached to the temples of the eyewear so it is out of the way.

The vertical angle of the light is adjustable so it can be pointed up or down depending where the object being viewed is located. To operate the device, press the on/off button once - this will turn on the full brightness of the LED. Press the button again to go to the lower level of light or Eco mode. Press it again and the device will turn off. A protective, hard foam case traveling case is included. The **Head Light Led** is available at both ViewFinder locations.



Visit us on the web at www.ViewFinderLowVision.com

1830 South Alma School Road, Suite 131 • Mesa, Arizona 85210 • 480.924.8755

10001 West Bell Road, Suite 115 • Sun City, Arizona 85351 • 623.583.2800

Our doctors also see patients at several satellite offices throughout the state of Arizona

Please call our toll free number 866.924.8755 or 480.924.8755

to schedule an appointment at one of these offices